




California State University of Bakersfield, Department of Chemistry

Ferrofluid



Standards:

3-PS2-3: Ask questions to determine cause and effect relationships of electric or magnetic interactions between two objects not in contact with each other.

3-PS2-4: Define a simple design problem that can be solved by applying scientific ideas about magnets.

MS-PS2-3: Ask questions about data to determine the factors that affect the strength of electric and magnetic forces.

MS-PS2-5: Conduct an investigation and evaluate the experimental design to provide evidence that fields exist between objects exerting forces on each other even though the objects are not in contact.

This material is based upon work supported by the Corporation for National and Community Service under Learn and Serve America Grant No. 10LHPCA001. Opinions or points of view expressed in this document are those of the authors and do not necessarily reflect the official position of the Corporation or the Learn and Serve America Program.

HS-PS3-5: Develop and use a model of two objects interacting through electric or magnetic fields to illustrate the forces between objects and the changes in energy of the objects due to the interaction.

Introduction:

This project involves making a magnetic liquid that is fun to play with and fun to make.

Materials:

- Iron III Chloride
- Iron II Chloride
- Ammonium Hydroxide
- Tetramethylammonium Hydroxide
- Hydrochloric Acid (1M)
- Mineral Oil
- Distilled Water
- Strong Magnet
- Stirring Rod
- 125 mL Erlenmeyer flask
- 50 mL Erlenmeyer flask
- Weighing Boat
- Small Pipets
- Liquid Waste Container
- Small Spatula
- Permanent (STRONG) Magnet

Safety:

- Always have an adult with you to help you during your experiment.
- Dispose any chemicals in a liquid waste container.
- Wear lab aprons, gloves, and goggles when performing this experiment.
- Hydrochloric acid can be highly corrosive. Handle with care.
- Keep all personal electronics away from work station.
- Be extremely careful with the magnet. Do not leave it around other metal objects.

Procedure:

Step 1

- Weigh 2.7g of Iron III salt using balance and weighing boat.
- Place iron salt into a 125 mL Erlenmeyer flask or beaker.
- Use graduated cylinder to measure 25 mL of distilled water.
- Use graduated cylinder to measure 10 mL of 1M hydrochloric acid.
- Add water, then HCL into iron and shake (swish) until the iron has been completely dissolved. **Note:** Iron may get stuck to bottom. Use a stirring rod to scrape off and mix.

Step 2

- Weigh 1g of Iron II salt as before with balance
- Place Iron salt in a small beaker container


This material is based upon work supported by the Corporation for National and Community Service under Learn and Serve America Grant No. 10LHPCA001. Opinions or points of view expressed in this document are those of the authors and do not necessarily reflect the official position of the Corporation or the Learn and Serve America Program.

- Measure 2.5 mL of 1M hydrochloric acid and add to the Iron salt
- Shake until iron salt is completely dissolved

Step 3

- Combine the iron solutions by adding the Iron II solution with the Iron III.
- Mix the solutions well by swishing.
- Wash out the Iron II beaker with a small amount of distilled water and add to the solution.
- Mix solution again by shaking
- Add 50 mL of ammonium hydroxide by drops while shaking as much as possible.
- After all is in the solution, shake for at least three minutes.

Step 4



- Put the magnet under the beaker
- Wait until upper solution gets clear
- Decant solution into a waste beaker.
- Take 50 mL distilled water and add to magnetic particles
- Repeat this washing step at least 2 more times

Step 5

- Leave 50 mL of distilled water in the solution and shake to make all black again.
- Transfer the solution in small amounts to a weighing boat, using the magnet to settle the particles.
- Decant the clear upper solution into the waste beaker. Hold the magnet when decanting. **Note:** To remove everything, extra water may be required.
- Once all magnetic particles are in the weighing boat, get rid of as much water as possible.

Step 6

- Measure 5 mL of 25% tetramethylammonium hydroxide by using a pipette and add it to the magnetic particles stir/mix with a glass rod. Wait at least 3 minutes. Remove the black liquid on top.

Step 7

This material is based upon work supported by the Corporation for National and Community Service under Learn and Serve America Grant No. 10LHPCA001. Opinions or points of view expressed in this document are those of the authors and do not necessarily reflect the official position of the Corporation or the Learn and Serve America Program.

- Check the ferrofluid with a magnet for spikes. If there are no spikes, there is too much water.
- Fill a vial with mineral oil and add just a few drops of ferrofluid.
- Seal the container; try to leave as much air out as possible.

Data and Observations:

Record your observations in this space. What did you see? Anything you were not expecting? Something really awesome? Describe it here.

References:

Make Your Own Ferrofluid. KoserLab. <https://sites.google.com/site/koserlab/outreach/make-your-own-ferrofluid> (Accessed July 18, 2013)